

## Triggered Assessment in Lung Ultrasound (Supervisor)

In each of the following 3 sections the learner must achieve	Assessors comments	Competent?	
		Yes	No
<b>1. Preparation</b>		Tick	Tick
Demonstration of appropriate attitude and professional manner		Tick	Tick
Explanation of procedure to patient as appropriate		Tick	Tick
Checking patient's details/entry into machine as appropriate		Tick	Tick
Confirmation of indication and check any supportive imaging		Tick	Tick
Positioning of patient and machine ergonomically		Tick	Tick
<b>2. The Scan</b>		Tick	Tick
Selection of appropriate probe(s) and optimisation of settings		Tick	Tick
Places probe at upper, lower and postero-lateral points bilaterally in a longitudinal (cephalad/caudad) plane		Tick	Tick
Performs a systematic examination		Tick	Tick
Identification of subcutaneous tissue, ribs and pleura		Tick	Tick
Identification of pleural sliding in 2D and M mode		Tick	Tick
Identification of A lines, B lines, diaphragm, liver and spleen		Tick	Tick
Identifies any abnormalities and pathology correctly		Tick	Tick
<b>3. Post scan</b>		Tick	Tick
Adequate documentation and storage of images as appropriate		Tick	Tick
Informing patient and reporting findings where appropriate		Tick	Tick
Identification whether need for further scan or alternative imaging		Tick	Tick
Cleaning of equipment and storage to minimise damage		Tick	Tick
<b>Overall performance, recording and interpretation of scan</b>		<b>Tick</b>	<b>Tick</b>